

# High Tibial Osteotomy System

**BIDDU**<sup>®</sup>

## Tibial A/P Sloped Wedge Osteotomy Plate

- Opening wedge.
- Precise, sloped plate.
- For tibial Osteotomy.
- Maintains normal anatomy.
- Available in sizes 7.5mm, 10mm, 12.5mm and 15mm.
- Material - Stainless Steel 316L.



### BIDDU<sup>®</sup> Plates and Screws (Size Options)

BAK-7094.075S	BIDDU <sup>®</sup> Tibial A / P Sloped Wedge Osteotomy Plate, 7.5mm
BAK-7094.100S	BIDDU <sup>®</sup> Tibial A / P Sloped Wedge Osteotomy Plate, 10.0mm
BAK-7094.125S	BIDDU <sup>®</sup> Tibial A / P Sloped Wedge Osteotomy Plate, 12.5mm
BAK-7094.150S	BIDDU <sup>®</sup> Tibial A / P Sloped Wedge Osteotomy Plate, 15.0mm
BAK-7181.45.28	BIDDU <sup>®</sup> LPH Cortex Screw, Self tapping, Dia. 4.5mm, Len. 28mm
BAK-7181.45.30	BIDDU <sup>®</sup> LPH Cortex Screw, self tapping, Dia. 4.5mm, Len. 30mm
BAK-7181.45.32	BIDDU <sup>®</sup> LPH Cortex Screw, self tapping, Dia. 4.5mm, Len. 32mm
BAK-7181.45.34	BIDDU <sup>®</sup> LPH Cortex Screw, self tapping, Dia. 4.5mm, Len. 34mm
BAK-7181.45.36	BIDDU <sup>®</sup> LPH Cortex Screw, self tapping, Dia. 4.5mm, Len. 36mm
BAK-7181.45.38	BIDDU <sup>®</sup> LPH Cortex Screw, self tapping, Dia. 4.5mm, Len. 38mm
BAK-7181.45.40	BIDDU <sup>®</sup> LPH Cortex Screw, self tapping, Dia. 4.5mm, Len. 40mm
BAK-7182.65.45	BIDDU <sup>®</sup> LPH Cancellous Screw, fully-threaded, Dia. 6.5mm, Len. 45mm
BAK-7182.65.50	BIDDU <sup>®</sup> LPH Cancellous Screw, fully-threaded, Dia. 6.5mm, Len. 50mm
BAK-7182.65.55	BIDDU <sup>®</sup> LPH Cancellous Screw, fully-threaded, Dia. 6.5mm, Len. 55mm
BAK-7182.65.60	BIDDU <sup>®</sup> LPH Cancellous Screw, fully-threaded, Dia. 6.5mm, Len. 60mm
BAK-7182.65.65	BIDDU <sup>®</sup> LPH Cancellous Screw, fully-threaded, Dia. 6.5mm, Len. 65mm
BAK-7182.65.70	BIDDU <sup>®</sup> LPH Cancellous Screw, fully-threaded, Dia. 6.5mm, Len. 70mm

## BIDDU<sup>®</sup> Instrument set



### Ordering Information

BITO-001	BIDDU <sup>®</sup> Bone Spreader
BITO-002	BIDDU <sup>®</sup> Osteotomy Chisel, 10mm
BITO-003	BIDDU <sup>®</sup> Osteotomy Chisel, 15mm
BITO-004	BIDDU <sup>®</sup> Osteotomy Chisel, 20mm
BITO-005	BIDDU <sup>®</sup> Osteotomy Chisel, 25mm
BITO-006	BIDDU <sup>®</sup> Plate Holder
BITO-007	BIDDU <sup>®</sup> Screw Driver SW 3,5 with Holding Sleeve
BITO-008	Drill bit Dia. 3,2mm; Len. 150mm (For BIDDU <sup>®</sup> LPH Cortex Screw, Dia. 4,5mm)
BITO-009	Drill bit Dia. 3,5mm, Len. 150mm (For BIDDU <sup>®</sup> LPH Cancellous Screw, Dia. 6,5mm)
BITO-010	BIDDU <sup>®</sup> Osteotomy Wedge Sizer
BITO-011	BIDDU <sup>®</sup> Depth Gauge, Large
BITO-012	BIDDU <sup>®</sup> Screw Holding Forceps
BITO-013	BIDDU <sup>®</sup> Osteotomy Jack
BITO-014	Graphic Case for BIDDU <sup>®</sup> Instrumentation Set
BITO-016	Alignment Rod, Male/Female
BITO-023	Gap Measuring Device

