

# Small Joints Ligament Repair

## **MICRO-VIM™** Ligament Anchor, Titanium

MICRO-VIM™ Ligament Anchor is made of radiopaque titanium. The screw-in design of anchor is ideal for soft-bone matrix. The MICRO-VIM™ Ligament Anchor is available in 1.8mm and 2.0mm diameter. The Inserter is preloaded with anchor and 1 pc USP size 2-0/0 Fiber with tapered needle. The slotted sleeve and anchor specific drill bits are available for accurate placement of the anchor.



### Ordering Information

BAS-9161.18F MICRO-VIM™ Ligament Anchor, Titanium, loaded with 1 pc Fiber USP 2-0 with needle, Dia. 1.8mm  
BAS-9162.20F MICRO-VIM™ Ligament Anchor, Titanium, loaded with 1 pc Fiber USP 0 with needle, Dia. 2.0 mm

### Related Instrumentation

BAS-9148 Drill bit, Dia. 1.5mm for MICRO-VIM™ Ligament Anchor BAS-9162.20F  
BAS-9149 Drill bit, Dia. 1.3mm for MICRO-VIM™ Ligament Anchor BAS-9161.18F  
BAS-9150 Slotted Sleeve for MICRO-VIM™ & MICRO-VIM™ PK Ligament Anchor  
BAS-9161.18F, BAS-9162.20F, BAS-9163.18F & BAS-9164.20F

# MICRO-VIM™ PK

## Ligament Anchor, PEEK

MICRO-VIM™ PK Ligament Anchor is made of radiolucent PEEK polymer. The tap-in design of anchor is ideal for hard-bone matrix. The MICRO-VIM™ PK Ligament Anchor is available in 1.8mm and 2.0mm diameter. The Inserter is preloaded with anchor and 1 pc USP size 2-0/0 Fiber with tapered needle. The slotted sleeve and anchor specific drill bits are available for accurate placement of the anchor.



### Ordering Information

- BAS-9163.18F MICRO-VIM™ PK Ligament Anchor, PEEK, loaded with 1 pc Fiber USP 2-0 with needle, Dia. 1.8mm
- BAS-9164.20F MICRO-VIM™ PK Ligament Anchor, PEEK, loaded with 1 pc Fiber USP 0 with needle, Dia. 2.0 mm

### Related Instrumentation

- BAS-9146 Drill bit, Dia. 1.7mm for MICRO-VIM™ PK Ligament Anchor BAS-9164.20F
- BAS-9147 Drill bit, Dia. 1.5mm for MICRO-VIM™ PK Ligament Anchor BAS-9163.18F
- BAS-9150 Slotted Sleeve for MICRO-VIM™ & MICRO-VIM™ PK Ligament Anchor  
BAS-9161.18F, BAS-9162.20F, BAS-9163.18F & BAS-9164.20F

## Indications for MICRO-VIM™ and MICRO-VIM™ PK Ligament Anchor

Hand/Wrist	Foot/Ankle
Scapholunate Ligament Reconstruction	Mid-foot reconstruction
TFCC repair	Metatarsal/Phalanges Ligament Repair
Repair/reconstruction of flexor and extensor tendon at the PIP (proximal interphalangeal), DIP (distal interphalangeal), and MCP (metacarpal interphalangeal) joints for all digits	Plantar Plate Repair
Ulnar or lateral collateral ligament reconstruction	

