

4.0mm Cannulated Cancellous Screw System

1.25 mm Non-Threaded Guide Wire, 150 mm



1.25 mm Threaded Guide Wire, 150 mm
Maintains reduction during drilling. Threaded spade point tip allows easy penetration into the bone and maximum resistance to inadvertent removal.



2.7 mm Cannulated Drill Bit, 160 mm, 1.35 mm
cannulation, quick coupling



Cannulated Hexagonal Screwdriver, 2.5 mm hex
Fully cannulated for insertion of any 4.0 mm
cannulated screw over the guide wire.
Accommodates holding sleeve.



Cannulated Screw Measuring Device Provides a direct reading for screw length. Measurement places screw tip at thread/shaft junction of guide wire.



1.25 mm Cleaning Stylet For intraoperative cleaning of 4.0 mm cannulated screw instruments.

