

PELVIC BASIC INSTRUMENT SET



The complex anatomy of the pelvic ring and acetabulum involve complex surgical exposures and demands exact anatomical reduction to achieve positive functional results.

Pelvic Standard Instrument Set

The Pelvic Instrument Set has a wide range of standard instruments for screw insertion, plate adaptation (bending instruments) and fixation during reduction for preliminary and definitive stabilization of pelvic and acetabular fractures.

Drill Bit, 3.5mm dia., L 195/170mm, for quick coupling



Tap for Cortex Screws Ø 3.5 mm, calibrated, length 175 mm



Holding Sleeve



Screwdriver, hexagonal, small, Ø 2.5 mm, length 270 mm



Drill Bit, 2.5mm dia., L 230/205mm, 3-flute, with depth marking



Drill Bit Ø 2.5 mm, calibrated, length 300 mm, with Quick Coupling, for Percutaneous Insertion



Depth Gauge for Cortex Screws Ø 3.5 mm, measuring range up to 150 mm



Bending Iron for Reconstruction Plates 3.5 and 4.5,
length 190 mm



Bending Pliers for Reconstruction Plates 3.5



Ratchet Wrench with AO/ASIF Quick Coupling
Handle with Quick Coupling



Screwdriver Shaft, hexagonal for
Screws \varnothing 3.5 mm, length 250 mm




Screwdriver Shaft **Stardrive**[®],
length 250mm



Screw Forceps, self-holding, length 85 mm



In-situ Bending and Twisting
Handle, straight



In-situ Bending and Twisting
Handle, 90°



In-situ Bending and Twisting
Handle, 120°



In-situ Bender




PELVIC REDUCTION & RETRACTOR INSTRUMENT SET

Pelvic Retractor Set

A wide range of retractors (pelvic, Hohmann, sciatic nerve and malleable) within the Pelvic Retractor Set provide more options to aid the surgeon in achieving good visibility within the chosen approach.

Pelvic Reduction Instrument Set

The pelvic reduction forceps and other reduction instruments within the Pelvic Reduction Instrument Set are designed to aid the surgeon with fracture reduction of pelvic injuries.



Schanz Screw Ø 6.0 mm, length 190/50 |
mm, Stainless Steel



Universal Chuck with T-Handle, L 160mm



Bone Hook, sharp, medium, L 230mm



Bone Hook, sharp, large, L 200mm



Pelvic Retractor, blunt, L 270mm



Pliers for Screw Removal, L 205mm



Pelvic Reduction Forceps, medium, length 250 mm, for use with Cortex Screws Ø 3.5 and 4.5 mm



Pelvic Reduction Forceps, large, adjustable, speed lock, length 330 mm



Retractor, 18mm wide, long narrow tip, L 235mm



Reduction Forceps, large, with points, L 200mm




Pelvic Retractor, medium, length 268 mm



Pelvic Retractor, large, length 323 mm




Retractor, can be contoured, width 20 mm



Retractor, can be contoured, width 30 mm



Retractor, can be contoured, width 40 mm



Ball Spike, straight, with pointed ball tips Ø 6.5 mm, length 300 mm



Ball Spike, straight, long, with pointed ball tips Ø 6.5 mm, length 400 mm



Pelvic Reduction Forceps with pointed ball tips Ø 6.5 mm, length 250 mm



Pelvic Reduction Forceps, oblique, with pointed ball tips Ø 6.5 mm, length 200 mm



Pelvic Reduction Forceps, oblique, with pointed ball tips \varnothing 6.5 mm, length 240 mm



Pelvic Reduction Forceps with pointed ball tips \varnothing 6.5 mm, length 400 mm



Pelvic Reduction Forceps, asymmetric, with pointed ball tips \varnothing 6.5 mm



Disc, hole \varnothing 6.5 mm



Disc, rectangular, hole \varnothing 6.5 mm



Bone Hook, medium, long, length 330 mm



Bone Hook, large, long, length 330 mm



Pelvic Reduction Forceps, small, length 190 mm, for use with Cortex Screws \varnothing 3.5 and 4.5 mm



Pelvic Reduction Forceps, small, for Screws \varnothing 3.5 mm, length 250 mm





BIOTEK - Chetan Meditech Pvt. Ltd.
MD-04, Sanand-II, GIDC,
Ahmedabad Gujarat INDIA- 382110
Email: contact@biotekortho.com
URL: biotekortho.com



An ISO 13485 : 2016 Company